



December 2020

xsens

Release notes MVN
2020.2

MVN 2020.2 release notes

New versions	MVN: 2020.2 (Dec 2020) Awinda firmware: (unchanged) AW-A2 v4.4.0, AW-DNG2 v4.4.0, MTw2 v4.4.0 Link firmware: BPACK v1.6.0, MTx2 v1.6.0, BPACKv2 v1.6.0 FWU: 4.3.1
Previous versions	MVN: 2020.0.2 (Sept 2020) Awinda firmware: AW-A2 v4.4.0, AW-DNG2 v4.4.0, MTw2 v4.4.0 Link firmware: BPACK v1.2.0, BPACKv2 v1.2.0 MTx2 v1.2.0 FWU: 4.2.7
Release date	December 2020
Operating system	Windows 10, 64-bit
Meaning of icons	+added, ^changed, !note, #fixed, -removed

1 Introduction

This document describes the new features, improvements and bug fixes in the new release of Xsens MVN Analyze/Animate 2020.2

NOTE: this release does require a firmware update of your MVN Link system.

2 Major improvements and new features

2.1 MVN

- + HTC Vive live position aiding added
- + HTC Vive live Object tracking added
- + support of BPv2 in MVN
- + Adjustable update rate for MVN Link
- + Network sync slave function
- + ROBR wireless download data
- + Data stream from StretchSense Glove
- + Live stream data in Matlab
- + EMG integration
- + Velocity and acceleration CoM
- + set body dimensions are mandatory at ROBR
- ^ Manus Prime II integration

2.2 Firmware

It is required to upgrade the firmware of the different components of your MVN Link to the latest versions to enjoy the latest features in MVN.

Download the Firmware Updater (FWU) from the Xsens website.

3 Bug fixes

3.1 MVN

- # Changing suit config changes mouse to a text cursor
- # Fixing uninstalls
- # Fix on Batch reprocessing, with empty MVN files
- # Exporter options are translated
- # Files creation no longer creates delay in the start of a recording
- # Checkboxes in splash screen now have same behaviour
- # Battery life of link now updated during recordings
- # MVN in Japanese can now find mvna files again
- # Video file now immediately synced
- # Fixed stepped data in Manus + MVN export
- # Fixed using CTRL + Z in contact editing

4 Minor updates

4.1 MVN

- BPv2 support
- Recording start/stop delays added for all licenses
- MVN can be used without license, recording only license
- Delsys EMG support
- MVN + EMG excel exporter added
- Smoothing filter add for object tracking
- Objects can be measured without connecting a MVN system
- ROBR UI updated
- StretchSense support
- EULA add to copy to clipboard in the About Xsens window
- Vive floor calibration
- Cometa EMG support
- Unity stream option renamed to Unity (legacy)
- HTC HMD can be used as position aiding
- HTC controller can be used as position aiding
- Up to 5 position aiding devices to 1 subject/actor
- Ergonomics joint angles in network streamer

5 Known Issues

-