

MVN Remote app 2.0 release notes

New versions	MVN Remote 2.0 Supported MVN version: MVN 2019.2.1
Previous versions	MVN Remote 1.3 Supported MVN version: MVN 4.3.9
Release date	November 2019
Meaning of icons	+added, ^changed, !note, #fixed, -removed

1 Introduction

This document describes the new MVN Remote app 2.0.
This release contains new features and improvements, as well as some bug fixes.

NOTE: To take advantage of the new features and stability updates, you need the supported MVN version (see above).

2 Major improvements and new features

- + Support new reference video handling in MVN Remote app
- + Add reference video transfer progress reporting
- + Add option to automatically delete reference videos from mobile device after transfer completes
- + Sync MVN software and MVN Remote clocks
- + Add timestamps to video recordings
- + Finger data support added in MVN Remote app
- + Prop visualization added in MVN Remote app

3 Bug fixes

- # Multiple person can be shown
- # Move character to Origin fixed
- # Preselected video resolution
- # Adjustable video resolution
- # Reconnection to MVN fix

4 Known Issues

- In Windows 7, the live video stream is not supported